

The logo for Integral Molecular, featuring the word "integral" in a large, bold, white sans-serif font, with "molecular" in a smaller, italicized white sans-serif font below it. The text is set against a dark blue background with a faint, stylized graphic of a protein structure or membrane protein.

**integral**  
molecular

Membrane Protein Solutions

# Harnessing Divergent Species to Access Difficult and Conserved Antibody Targets


---

*April 25, 2024*  
*Antibody Society Webinar*

Ross Chambers PhD  
VP of Antibody Discovery  
Integral Molecular, Philadelphia PA

# Integral Molecular

The Industry Leader in Delivering Lead Antibodies Against Undruggable Targets





Delivering MAbs against:

- Complex membrane proteins
- Highly conserved epitopes
- Functional epitopes
- Epitopes with rare properties

- 20+ years working with challenging protein targets
- Industry-leading epitope diversity
- Pipeline of therapeutic antibodies against complex targets
- Therapeutic MAbs licensed to AstraZeneca, Context Therapeutics, Cartexell, and others

# Trusted by 600+ Companies

500+ publications and patents, including *Cell*, *Science*, *Nature*



# Integral Molecular

The Industry Leader in Delivering Lead Antibodies Against Undruggable Targets




## Breaking barriers in antibody discovery:

harnessing divergent species for accessing difficult and conserved targets


Banik, S.R., Kushnir, N., Doranz, B.J., and Chambers, R.

### ABSTRACT

To exploit highly conserved and difficult targets, including multipass membrane proteins, monoclonal antibody discovery efforts increasingly rely on the advantages offered by divergent species such as rabbits, camelids, and chickens. Here, we provide an overview of antibody discovery technologies, analyze gaps in therapeutic antibodies that stem from the historic use of mice, and examine opportunities to exploit previously inaccessible targets through discovery now possible in alternate species. We summarize the clinical development of antibodies raised from divergent species, discussing how these animals enable robust immune responses against highly conserved binding sites and enable antibodies capable of penetrating functional pockets via long HCDR3 regions. We also discuss the value of pan-reactive molecules often produced by these hosts, and how these antibodies can be tested in accessible animal models, offering a faster path to clinical development.



Conservation of Drug Targets in Chickens



December 2023

# How Divergent Species Can Access Conserved Targets



Gaps in antibody space & role of divergent species



Rabbits, Camelids, Chickens




Chicken immunization has delivered:


Antibodies against conserved targets  
Agonist antibody  
Exquisite specificity

# Why Are Conserved Targets Difficult for Antibody Discovery?

## Robust Immune Reaction




## Immune Tolerance




# Why Are Conserved Targets Difficult for Antibody Discovery?

## Robust Immune Reaction




## Immune Tolerance




# Why Are Conserved Targets Difficult for Antibody Discovery?

## Robust Immune Reaction



## Immune Tolerance




- Immune system recognizes self-antigen
  - Tolerance prevents immune system from attacking own body
- Conserved proteins resemble self-antigen

Conserved proteins make poor immunogens for therapeutic discovery




# Many Highly Conserved Targets Currently Undruggable by MAbs



Banik et al., 2023, *mAbs*

Untapped opportunities for MAb discovery


# ~Half of Drug Targets May Be Poor Immunogens in Mice



# Antibody Enabling Technologies and Therapeutics


## Enabling Technologies

## Therapeutic Antibody Development



# Divergent Species Increasingly Used for MAb Discovery


## MAbs From Divergent Species in Clinical and Preclinical Development





## Advantages of Divergent Species

1. Less immune tolerance
2. Long HCDR3 regions
3. Cross species-reactive antibodies


# 1. Divergent Species Avoid Immune Tolerance



# 1. Divergent Species Avoid Immune Tolerance



# 1. Divergent Species Avoid Immune Tolerance



- Target antigen must be 'foreign' for host to mount effective immune response
- Sequence divergence important at epitopes of interest
- Chickens: long evolutionary distance + MAb structure similar to humans


Divergent species enable access to more targets and more epitopes

## 2. Longer HCDR3 Can Form Protruding Structures

- HCDR3 region correlates with paratope shape
- Ramsland et al., studied 50 antibody crystal structures


### HCDR3:

- **1-9 aa:** predominantly non-protruding
- **10-13 aa:** variable
- **14+ aa:** predominantly protruding






# Longer HCDR3 Can Access Pockets



# Divergent Species Can Access Recessed Epitopes


- HCDR3 of antibody is most important for contacting the epitope
- Mouse HCDR3 relatively short
  - Tend to bind flat epitopes
- Longer sequences can protrude into functionally important recessed epitopes



*\*Kabat numbering used throughout*

### 3. Cross Species-Reactivity Could Expedite MAb Development

- Many well-validated preclinical animal models
- Most mouse-derived MAbs reactive only in primates
  - Testing limited to non-human primates (NHP)
  - NHP experiments resource intensive
  - NHPs in short supply
  - NHP testing not actually required
- Ability to test MAbs in lower mammals could greatly streamline studies




#### Mammalian species used in animal models

\*in Nobel prize work (physiology/medicine)  
Adapted from (Jota Baptista et al., 2021, Pharmacology)

### 3. Cross Species-Reactivity Could Expedite MAb Development

- Many well-validated preclinical animal models
- Most mouse-derived MAbs reactive only in primates
  - Testing limited to non-human primates (NHP)
  - NHP experiments resource intensive
  - NHPs in short supply
  - NHP testing not actually required
- Ability to test MAbs in lower mammals could greatly streamline studies



#### Mammalian species used in animal models

\*in Nobel prize work (physiology/medicine)  
Adapted from (Jota Baptista et al., 2021, Pharmacology)

# Cross Species-Reactive MAbs From Divergent Species


	Clinical Stage	Immunization host	% Identity vs host	% Identity vs mouse	Reactivity in non-primates
<b>Caplacizumab</b> (von Willebrand factor)	Approved	Camelid	81	83	Guinea pig, mini-pig, pig,
<b>Brolucizumab</b> (VEGF-A)	Approved	Rabbit	88	84	Rat, mouse, dog, pig, cat
<b>Eptinezumab</b> (CGRP)	Approved	Rabbit	89	89	Rat, rabbit
<b>Sym021</b> (PD1)	Clinical trials	Chicken	35	59	Mouse
<b>CTIM-76</b> (Claudin 6)	Preclinical	Chicken	61	88	Mouse
<b>IM-68-27G10</b> (GPRC5D)	Preclinical	Chicken	49	82	Mouse
<b>B30</b> (BDNF)	Discovery	Chicken	94	100	Rat
<b>pT231/pS235_1</b> (Tau)	Discovery	Chicken	84	90	Mouse
<b>AC1</b> (CD20)	Discovery	Chicken	No ortholog	75	Mouse
<b>YW33</b> (Integrin $\alpha 11\beta 1$ )	Discovery	Chicken	78/85	90/93	Mouse, rat
<b>MAb panel</b> (SIRP $\alpha$ )	Discovery	Chicken	42	66	Mouse
<b>MAb panel</b> (SLC2A4)	Discovery	Chicken	65 (Paralog)	95	Mouse
<b>MAb panel</b> (GIPR)	Discovery	Chicken	50	82	Mouse, rat
<b>MAb panel</b> (Kv1.3)	Discovery	Chicken	84	96	Mouse
<b>MAb panel</b> (Claudin18.2)	Discovery	Chicken	76	90	Mouse

Banik et al., 2023, *mAbs*

# Our Experience Generating Species Cross-Reactive MAbs


## CLDN18.2

(90% identical to mouse)




## GPRC5D

(82% identical to mouse)



## CCR8

(71% identical to mouse)



(Antibodies raised in chickens)

# How Divergent Species Can Access Conserved Targets



Gaps in antibody space &  
role of divergent species



**Rabbits, Camelids, Chickens**



**Chicken immunization has delivered:**

Antibodies against conserved targets  
Agonist antibody  
Exquisite specificity

# Rabbit MAb Discovery

- Well-established role as immunological hosts
- Single VH (VH1) and VL (VK1) framework
  - Cloning, humanization, and engineering relatively straightforward
- MAb discovery technologies
  - B cell cloning (e.g. eptinezumab, crovalimab, clazakizumab)
  - Rabbit hybridomas (e.g. APX005M/sotigalimab)

	Mouse	Rabbit
Evolutionary distance from humans, years	91 million	91 million
Animal host class	mammalia	mammalia
Robust immune response against conserved proteins	-	-
Cross-reactive MAbs for preclinical models	-	+/-
Canonical IgG	+	+
HCDR3 length (aa)	9.1±2.0	12.8±3.6
Long paratope with average HCDR3 >14 aa	-	-
Nanobodies from host	-	-
Immunization costs and animal logistics	\$	\$
Access to humanized animals	+	-
Diverse accessible B cell repertoire (spleen, marrow)	+	+




# Rabbit MAb

INN	Brand name	Other names	Most advanced phase	Target(s)	Indications of clinical trials
Eptinezumab	Vyepti	ALD403	Approved in US and EU	Calcitonin gene-related peptide (CGRP)	Migraine prevention
Brolucizumab	Beovu	RTH258, ESBA1008, DLX-1008	Approved in US and EU	VEGF-A	Diabetic macular edema, neovascular age-related macular degeneration
Crovalimab		SKY59, RG6107, RO7112689	Regulatory review in US, EU, China and Japan	Complement C5	Sickle cell disease, atypical hemolytic uremic syndrome, paroxysmal nocturnal hemoglobinuria
Suvemcitug		APX003, BD0801; TK001, Sevacizumab (not WHO assigned INN)	Phase 3	VEGF	Ovarian cancer, macular degeneration
Clazakizumab		CSL300, ALD518; BMS-945429	Phase 3	IL-6	Subjects with end stage kidney disease undergoing dialysis, COVID-19, acute GvHD, psoriatic arthritis; Crohn's disease, rheumatoid arthritis, oral mucositis, non-small cell lung cancer-related fatigue and cachexia
		9MW0211	Phase 2/3	VEGF	Macular degeneration
Sotigalimab		APX005M	Phase 2 (non-commercial sponsor)	CD40	Ovarian cancer, melanoma, gastro-esophageal cancer, pancreatic cancer, solid tumors, non-small cell lung cancer, renal cancer
		LU AG09222, ALD1910	Phase 2	PACAP-38	Allergic rhinitis, migraine
		TRK-950	Phase 2 pending	CAPRIN-1	Gastric cancer, solid tumors
		YYB101	Phase 1/2	HGF	Colorectal cancer
		ASKB589	Phase 1/2	Claudin 18.2	Solid tumors
		CLM-101, NOV-110501, YYB-101	Phase 1/2	HGF	Colorectal cancer, solid tumors
QX005N, SNC005	Phase 1	IL-4R	Atopic dermatitis		

(Antibody Society, December 2023)

# Camelid MAb Discovery

Canonical MAb



VHH nanobody



- Small heavy-chain only MAbs
  - Can be engineered into new formats (bispecifics)
  - Long HCDR3
- MAb discovery techniques
  - Phage display (caplacizumab, envafolimab, and ozoralizumab)
  - Yeast display
- Complex logistics due to large host animals

	Mouse	Camelid
Evolutionary distance from humans, years	91 million	92 million
Animal host class	mammalia	mammalia
Robust immune response against conserved proteins	-	-
Cross-reactive MAbs for preclinical models	-	+/-
Canonical IgG	+	+
HCDR3 length (aa)	9.1±2.0	14.5±3.7
Long paratope with average HCDR3 >14 aa	-	+
Nanobodies from host	-	+
Immunization costs and animal logistics	\$	\$\$\$
Access to humanized animals	+	-
Diverse accessible B cell repertoire (spleen, marrow)	+	-


# Camelid MAbs

INN	Brand name	Other names	Most advanced phase	Target(s)	Indications of clinical trials
Caplacizumab	Cablivi	ALX-0081 (IV), ALX-0681	Approved in EU and US	von Willebrand factor	Acquired thrombotic thrombocytopenic purpura, unstable angina; non ST segment elevation myocardial infarction (NSTEMI); stable angina (associated with high risk PCI), thrombotic thrombocytopenic purpura
Ozoralizumab		TS-152, ATN-103	Approved in Japan	TNF, albumin	Rheumatoid arthritis
Envafoлимab	ENWEIDA	KN035, ASC22	Approved in China	PD-L1	HIV infection, metastatic or recurrent non-microsatellite highly unstable (non-MSI-H)/non-DNA mismatch repair defect (non-dMMR) endometrial cancer, undifferentiated pleomorphic sarcoma, hepatitis B, bile tract carcinoma, solid tumors
Netakimab	Efleira	BCD-085	Approved in Russia	IL-17	Ankylosing spondylitis, psoriasis, psoriatic arthritis
Levilimab	Ilsira	BCD-089	Approved in Russia	IL-6R	COVID-19, rheumatoid arthritis
Gefurulumab		ALXN1720	Phase 3	C5	Myasthenia gravis
Erfonrilimab		KN046	Phase 3	PD-L1/CTLA-4	Multiple tumor types
		PM8002	Phase 2/3	PD-L1, VEGF	Hepatocellular carcinoma, small cell lung cancer, non-small cell lung cancer, solid tumors
		LMN-201	Phase 2/3 pending	C. difficile exotoxin B	Clostridioides difficile infection
Ozekibart		JCT205, INBRX-109	Phase 2, pivotal	DR5	Chondrosarcoma, solid tumors including sarcomas
		BI 836880	Phase 2	VEGF, Ang2	Head and neck cancer, liver cancer, anal canal squamous cell carcinoma, macular degeneration, non-small cell lung cancer, solid tumors
		PM8001	Phase 2	PD-L1, TGFβ	Cancer
Sonelokimab		M1095, MSB0010841, ALX-0761	Phase 2	IL-17A, IL17F, HSA	Psoriatic arthritis, hidradenitis suppurativa, psoriasis
		LMN-101	Phase 2	FlaA	Campylobacter jejuni infection
Livmoniplimab		ABBV-151, ARGX-115	Phase 2	GARP-TGFβ1 complex	Solid tumors
Cusatuzumab		ARGX-110, JNJ-74494550	Phase 2	CD70	Cutaneous T-cell lymphoma, AML, hematological and solid cancers, Waldenström's macroglobulinemia
		LEO 138559, ARGX112, LP0145	Phase 2	IL-22R	Atopic dermatitis
		SAR442970	Phase 2	TNF, OX40L	Hidradenitis suppurativa
Tarperprumig		ALXN1820	Phase 2	Properdin	Sickle cell disease

+ many more early stage  
(Antibody Society, 2023)

# Chicken MAb Discovery

- Many conserved targets immunologically accessible
  - Only 15% of drug targets appear highly conserved (>90% identity)
- Canonical antibodies despite phylogenetic distance
- Skewed distribution of HCDR3s with longer sequences
  - ~90% with HCDR3  $\geq 13$ + amino acids
- Easy to humanize, only 1 germline gene for heavy and light chain
- MAb discovery techniques:
  - Phage display (numerous preclinical examples)
  - B cell cloning (Sym021)
- First chicken MAb in clinic
  - Sym021, targeting PD-1 (Cyno and mouse cross reactive)



# Chicken MAb Discovery

- Many conserved targets immunologically accessible
  - Only 15% of drug targets appear highly conserved (>90% identity)
- Canonical antibodies despite phylogenetic distance
- Skewed distribution of HCDR3s with longer sequences
  - ~90% with HCDR3  $\geq 13$ + amino acids
- Easy to humanize, only 1 germline gene for heavy and light chain
- MAb discovery techniques:
  - Phage display (numerous preclinical examples)
  - B cell cloning (Sym021)
- First chicken MAb in clinic
  - Sym021, targeting PD-1 (Cyno and mouse cross reactive)

	Mouse	Chicken
<b>Evolutionary distance from humans, years</b>	91 million	310 million
<b>Animal host class</b>	mammalia	aves
<b>Robust immune response against conserved proteins</b>	-	+
<b>Cross-reactive MAbs for preclinical models</b>	-	+
<b>Canonical IgG</b>	+	+
<b>HCDR3 length (aa)</b>	9.1 $\pm$ 2.0	15.9 $\pm$ 2.8
<b>Long paratope with average HCDR3 &gt;14 aa</b>	-	+
<b>Nanobodies from host</b>	-	-
<b>Immunization costs and animal logistics</b>	\$	\$
<b>Access to humanized animals</b>	+	+
<b>Diverse accessible B cell repertoire (spleen, marrow)</b>	+	+

# How Divergent Species Can Access Conserved Targets



Gaps in antibody space &  
role of divergent species



Rabbits, Camelids, Chickens




**Chicken immunization has delivered:**

Antibodies against conserved targets  
Agonist antibody  
Exquisite specificity

# Chickens Are Integral to Our Platform

## Native Antigen




- Pioneers in DNA & mRNA immunization for MAb discovery
- Inventors of Lipoparticle technology
- 20+ years membrane protein expertise

## Chicken MAb Discovery




- Access conserved & difficult epitopes
- Diverse antibody candidate panels
- Humanized, pM affinity, developable

## Preclinical Leads







- Bispecific/therapeutic format
- Functional POC *in vitro* & *in vivo*
- Deliverable: 12-18 months to IND

# High-Titer Responses to Conserved Targets



	CLDN6	Kv1.3	CB1	SLC2A4
% Mouse identity	88	96	97	95
% Chicken identity	61	84	94	65 (paralog)

Serum titer



Chickens deliver diverse MAbs against wide range of targets, with **95% success**






# Claudin 6 Program for Oncology – Enabled by Chickens

## Challenges

Only 3 extracellular residues different from CLDN9 (widely expressed)

Structurally complex antigen



ECL identity:  
95%

Total aa identity:  
88%



## Therapeutic Potential

- Expressed in ovarian, NSC lung, teratomas, gastric tumors
- Not expressed in normal adult tissues





## Current Status

- Configured as a bispecific
- Potent antitumor effects in animals, and good developability
- IND filed, entering clinical trials mid-2024

# Antibody Diversity From Chicken Immunization








# CLND6 MAb Discovery: A Search for Specificity




- MPA encompasses human membrane proteome including 24 claudin family members
- Specificity profiling tool under review by FDA, for qualification as a Drug Development Tool

# CLND6 MAb Discovery: A Search for Specificity




Scrceni et al., 2022, *iScience*

# GPCR Agonist MAb Enabled by Chickens



- Therapeutic antibody discovery project completed for Merck
- Merck sought an agonist antibody



- ✓ Human
- ✓ Non-human primate
- ✓ Mouse
- ✓ Rat
- ✓ Dog

“ The unique CDR diversity generated in chickens provided an interaction capable of activating a GPCR that no other rodent or human phage display derived antibodies could do even to the same epitope region. ”

-Director of Antibody Discovery, Merck



# SLC2A4 MAb: Unique Activity and Long HCDR3

**PNAS**  
Proceedings of the National Academy of Sciences of the United States of America www.pnas.org

**Isolation of state-dependent monoclonal antibodies against the 12-transmembrane domain glucose transporter 4 using virus-like particles**

David F. Tucker, Jonathan T. Sullivan, Kimberly-Anne Mattia, Christine R. Fisher, Trevor Barnes, Manu N. Mabilla, Rona Wif, Chidananda Sulli, Meghan Pitts, Riley J. Payne, Moniquetta Hall, Duncan Huston-Paterson, Xiaoxiang Deng, Edgar Davidson, Sharon H. Willis, Benjamin J. Doranz, Ross Chambers, Joseph B. Rucker  
Integral Molecular, Philadelphia, PA


**Abstract**  
The insulin-responsive 12 transmembrane transporter GLUT4 changes conformation between an inward-open state and an outward-open state to actively facilitate cellular glucose uptake. Because of the difficulties of generating conformational MAbs against complex and highly-conserved membrane proteins, no reliable tools exist to measure GLUT4 at the cell surface, follow its trafficking, or detect the conformational state of the protein. Here we report the isolation and characterization of novel conformational monoclonal antibodies (MAbs) that recognize the extracellular and intracellular domains of GLUT4, including MAbs that are specific for the inward-open and outward-open states of GLUT4. MAbs against GLUT4 were generated using virus-like particles (VLPs) to present this complex membrane protein in its native conformation, and using a divergent host species (chickens) for immunization to overcome immune tolerance. As a result, the isolated MAbs recognize conformational epitopes on native GLUT4 in cells, with apparent affinities as high as 1 pM and with specificity for GLUT4 across the human membrane proteome. Epitope mapping using shotgun mutagenesis alanine scanning across the 509 amino acids of GLUT4 identified the binding epitopes for MAbs specific for the states of GLUT4, as well as comprehensive identification of the residues that functionally control the GLUT4 inward-open and outward-open states. The MAbs identified here will be valuable molecular tools for monitoring GLUT4 structure, function, and trafficking, for differentiating GLUT4 conformational states, and for the development of novel therapeutics for the treatment of diabetes.




## SLC2A4 State-Specific MAb (LM048)

- 26 aa HCDR3
- Binds in cavity formed by ECL4 and ECL6
- Rare property of binding state-specific conformation

# Chicken MAbs Have Long HCDR3






# Chicken MAbs Have Long HCDR3





# Integral Molecular: Preclinical Pipeline of Chicken-Derived MAbs

Target	Indication	Modality	Discovery	Preclinical	Phase 1	Partner
<b>CLDN6</b>	solid tumors	TCE	Partnered			
<b>Undisclosed</b>	oncology	undisclosed	Partnered			
<b>CLDN18.2</b>	solid tumors	CAR-T	Partnered			
<b>CLDN18.2</b>	solid tumors	ADC	Internal (available)			Undisclosed
<b>CLDN18.2</b>	solid tumors	mRNA TCE	Internal (available)			
<b>GPRC5D</b>	multiple myeloma	TCE	Internal (available)			
<b>CCR8</b>	solid tumors	ADC	Internal (available)			
<b>KV1.3</b>	autoimmune	multiple	Internal (available)			

TCE: T- cell engager (multi-specific antibody)  
 ADC: Antibody-drug conjugate  
 CAR-T: Chimeric antigen receptor T cell



# Divergent Species Enable MAb Discovery



- Many valuable drug targets unfeasible in mice due to conservation
- Divergent species are increasingly used

# Divergent Species Enable MAb Discovery



- Many valuable drug targets unfeasible in mice due to conservation
- Divergent species are increasingly used

Evolutionarily divergent species provide:



- Robust immune response with broad epitope coverage
- Long HCDR3 regions can access functional pockets
- Access to more animal models because of species cross reactivity
- See Banik et al., 2023, mAbs

# Divergent Species Enable MAb Discovery



- Many valuable drug targets unfeasible in mice due to conservation
- Divergent species are increasingly used

Evolutionarily divergent species enable



- Robust immune response with broad epitope coverage
- Long HCDR3 regions can access functional pockets
- Access to more animal models because of species cross reactivity
- See Banik et al., 2023, *mAbs*



Chicken MAbs have demonstrated

- Exquisite epitope specificity
- Agonist & state specific activity

# Access Our Antibodies & Technologies



## Work With Integral Molecular

- Specificity Screening, Lipoparticles, Epitope Mapping (fee-for-service)
- Antibody discovery (partnerships)
- Therapeutic antibodies (licenses)



## VeRSaMAb Antibodies From Cell Surface Bio

- Validated Recombinant antibodies with unparalleled Specificity

*Tell us what M Ab you wish you had!*



**integral**  
*molecular*

Membrane Protein Solutions



## THANK YOU

Ross Chambers, PhD

VP of Antibody Discovery

[rchambers@integralmolecular.com](mailto:rchambers@integralmolecular.com)

[www.integralmolecular.com](http://www.integralmolecular.com)